

# Meetbuisjes

## Tenax Sorbent Tubes

Sorbent (coating)	Size (mm) ODxLength	Sections	Sorbent (mg)	Ends	Tube Cover	Separators	Pack Size	Part Number
Anasorb 747/Tenax TA	16x125	2	5.2 gm/1.2 gm	GO	-	WW	Each	226-171 ‡
PUF/Tenax/PUF	22x100	3	3cm/750mg/3cm	GO	P	WWW	Each	226-124*
Tenax	8x110	2	100/50	GS	B	WWW	10	226-35031
Tenax	8x110	2	100/50	GS	B	WWW	50	226-35-03
Tenax	6x70	2	20-okt	GO	A		50	226-35-01 *
Tenax	6x70	2	30/15	GS	A	FWW	50	226-35
Tenax GR/Anasorb GCB1	1/4x3-1/2 in	2	175/150	SS	B	SWS	Each	226-348 ‡ §
Tenax GR/Anasorb GCB1	1/4x3-1/2 in	2	125/120	GO	B	WWW	Each	226-345 ‡ §
Tenax TA	1/4x3-1/2 in	1	250	GO	B	WW	Each	226-360 ‡ §
Tenax TA	1/4x3-1/2 in	1	250	SS	B	SWWS	Each	226-357 ‡ §
Tenax TA	1/4x3-1/2 in	1	100	SS	B	SSW	Each	226-340 ‡
Tenax TA	1/4x3-1/2 in	1	100	GO	B	WW	Each	226-339 ‡
Tenax TA	16x125	1	1.6 gm	GO	-	WW	Each	226-134 *
Tenax/Glass Fiber Filter (OVS)	13 to 8x75	2	140/70	GO	V	FFGT	10	226-56 *
Tenax	6x130	1	17	GO	-	WW	50 sets	226-35-02 *

### Table Notes:

Catalog Number:

\* : Limited shelf-life; contact SKC for more information

‡ : Tubes are chemically conditioned before shipping; use within 30 days or recondition. Restocking fee applies.

§ : Each tube has a flow direction arrow and unique number.

• : Limited shelf-life; refrigerator/freezer storage may be required. Contact SKC for more information.

Tube Ends:

GS: Glass Sealed

GO: Glass Open

SS: Stainless Steel Open

Protective Tube Covers:

A: 6-mm OD x 70-mm length

B: 8-mm OD x 110-mm length

C: 10-mm OD x 150-mm length

D: 10-mm OD x 220-mm length and shorter tubes

P: Low-volume PUF Tube Holder Do not use an adjustable low flow holder

V: OVS Tube Holder Do not use an adjustable low flow holder

Meetapparatuur voor professionals



Separators:

W: Glass Wool

G: Glass Fiber Filter

F: Foam

N: Nylon Ring

Q: Quartz Filter

R: Glass Spacer

S: Screen

T: PTFE Ring

Sorbent Equivalencies:

Anasorb 708 is equivalent to Chromosorb 108.

Anasorb 727 is equivalent to Chromosorb 106.

Anasorb 747 is an SKC proprietary sorbent.

Anasorb C300 is equivalent to Hydrar, Carulite.

Anasorb CSC is an SKC proprietary sorbent.

Anasorb GCB1 is equivalent to Carbotrap B (20/40), Carbopack B (60/80).

Trademarks:

Anasorb ® — SKC Inc

Chromosorb ® — Manville Corp.

Porapak ® — Waters Associates, Inc.

Tenax ® — Buchem B.V.

XAD ® — Rohm & Haas

Meetapparatuur voor professionals

**ProCare** | SAFETY  
*We take intensive Care*

[www.procaresafety.nl](http://www.procaresafety.nl) | [info@procaresafety.nl](mailto:info@procaresafety.nl)

**Locatie Groningen**

Zernikepark 16a  
9747 AN Groningen  
T 050 2072172

**Locatie Amsterdam**

Olympisch Stadion 24-2  
1076 DE Amsterdam  
T 020 2613999

 @ProCare\_safety

